IMPORTANT: Please refer to the applicable current standard (AS/NZ/ISO) your practice is working under for the testing of reusable instruments, as well as the instructions for use, safety, maintenance & testing set out by the manufacturer of your specific ultrasonic cleaner (UC).



This kit is manufactured to AS2773.2-1999 for Ultrasonic Cleaners for health care facilities, Benchtop for the 6.2 (b)(ii) Pencil Load.

Inside you will find directions for use and a table that could be used to record your testing.

REORDER CODE: 9S-HEX www.stshealth.com.au

## **STS HEX TEST**

## Pencil Load Kit for Ultrasonic Cleaners to AS2773.2



REORDER CODE: 9S-HEX www.stshealth.com.au

DATE	PASS/FAIL	NAME	INITIAL	DATE	PASS/FAIL	PASS/FAIL NAME
					_	
					_	

IMPORTANT: Always follow the directions for use and safety set out by the manufacturer of the ultrasonic cleaner. Do not touch the HEX while UC is operating. Be careful to check water temperature at end of 3 minutes. After each use, please check the HEX surface and replace if smooth or when approaching 60 tests.

DIRECTIONS FOR USE AS2773.2-1999 6.2(b): Start with a clean tank, add water to operating level & detergent, degas, rub surface of HEX with supplied HB pencil lead, immerse disc on rack. The HEX should be completely free from lead within 3 minutes. When complete, remove HEX, record results, dry off & replace in bag with pencil. Before use, the key ring can be added to disk to assist with positioning in tank.