

## INCISION & COAGULATION

30 watts is enough for greater efficiency.

### Technological performances

- Controlled power
- Reliable and reproducible setting
- Safer with maximal efficacy
- Less energy loss

Controlled power for all electrosurgical applications

- Incisions/Excisions
- Frenectomy
- Gingivoplasty
- Coagulation
- Abscess incision
- Exposure of impacted/retained teeth
- Gingivectomy
- Etc.

### Conductive bracelet



ACTEON®'s choice for a bracelet system is the best compromise for safety and convenience:

- Easy to install to enhance patient acceptance
- In direct contact with the patient's skin
- Easily disinfected
- Adjusted to a contact region with a low resistance
- Perfect size to avoid concentration of heat
- More security: in compliance with new standards

### Reliable and compact unit

- Very compact
- User-friendly for a fast setting up
- Only two potentiometers (incision and coagulation)
- Settings can be very finely adjusted according to the different types of tissues encountered

### Electrode holder

- Light and autoclavable
- Simple-to-mount electrodes
- The electrode is insulated by a sheath, only the active part is in contact with the soft tissue




### Classification of electrodes according to their size and hemostatic capacity

Incision  
Ø 0.22mm



Excision  
Ø 0.22mm



Coagulating incision  
Ø 0.40mm




Fulguration and coagulation  
Ø 1mm / Ø 2.5mm / Ø 3.2mm




# TECHNIQUES OF USE

# CLINICAL CASES




## Electrosection

Vaporizing the tissue with the active electrode makes it possible to cut quickly without damaging the edges of the incision.





## Fulguration

Makes it possible to surface coagulate gingival tissue while protecting the underlying layer.




## Coagulation

Concentrating energy at the surface of a massive electrode makes it possible to diffuse the heat to the surrounding tissue and produce instantaneous hemostasis.



## Gingival Eviction

One of Servotome®'s best recommended uses; can easily be combined with the Expsasl® technique (Acteon® Pharma).




## Gingivectomy

For extremely precise tracing of the ideal gingival line.



## Abscess incision and drainage

Fast and deep treatment.



## Exposure of impacted tooth

Circular incision for a clean and successful exposure.

	SURGICAL APPLICATION	ELECTRODE(S) RECOMMENDED									
		I22S	I22CA	TR22T	TR22R	TR22L	I40S	I40CA	FC10N	FC25B	FC32B
ORAL SURGERY	Gingival eviction	■	■								
	Abscess incision and drainage	■	■				■	■			
	Exposure of impacted tooth - operculectomy	■	■								
	Removal of the pericoronal sac	■	■								
	Hemostasis						■	■	■	■	■
RESTORATIVE DENTISTRY	Dental neck exposure	■	■								
	Exposure of fractured root crown lengthening	■	■				■	■			
	Elimination of hypertrophic gingiva			□	□	□	■	■			
	Widening of sulcus	■	■				■	■	■		
ORTHODONTICS	Gingival modeling of edentulous crest			□			■	■			
	Exposure of impacted tooth	■	■								
	Frenectomy	■	■				■	■			
OPERATORY DENTISTRY	Crown lengthening for brackets placement	■	■				■	■			
	Elimination of hypertrophic gingiva	■	■	□	□	□	■	■		■	
	Frenectomy	■	■				■	■			
PERIODONTICS	Preparation of cavity for inlay before impression						■	■	■		
	Flap surgery	■	■								
	Gingivectomy	■	■	□	□		■	■			
	Gingivoplasty	■	■	□	□		■	■			
	Stripping	■	■	□	□		■	■			