

## **DIAO.** Outstanding performance – excellent control.

For more efficiency in everyday practice.



2 | **DIAO** 

## DIAO. The rose-gold standard in crown, bridge and cavity preparation.



#### Sharpness means control:

The outstanding initial cutting performance ensures optimum control and superb preparation results.

Concentrated power for exceptional durability Extremely long lasting sharpness



**Facilitates the daily practice routine:** Just reach out for brilliant rose gold

**DIAO** | 3

On average,

34% longer service life, an 27%

greater initial sharpness compared to traditional instruments\*

\*Source: Komet Dental TestLab, mechanical cutting test 2020

### Efficiency, created from diamonds and pearls.

/unl

The innovative diamond coating of the DIAO burs is interspersed with ceramic pearls. Like this, the pressure exerted is concentrated on the tips of the diamond grains which greatly increases their efficiency.









# **DIAO.** More efficiency in everyday practice.

Users rave about the innovative DIAO coating for crown, bridge and cavity preparation. DIAO instruments are easy to control and retain their sharpness over a very long period of time, making treatment more effective and less tiring, and the daily workflow at the practice more efficient.

The unique DIAO rose-gold shade makes the instrument easily recognizable so it can be quickly accessed for the next treatment, making the workflow at the practice by far more efficient. What makes the innovation of DIAO so unique? The combination of ceramic pearls and diamond grains ensures unparalleled performance and control of the instrument. The innovative diamond coating of DIAO provides a significantly longer service life than comparable instruments and also makes them easier to clean.



Comparison of durability



Komet

## **DIAO.** Innovation you can feel.

#### Concentrated cutting performance for optimum control and long lasting, effective substance removal

DIAO is much easier to handle and impresses the user with its smooth operation and the way it glides over the tooth. Although easy to control, DIAO instruments remove substance consistently and at greater speed.







#### Directions: for use:

 Use at an optimum speed of O<sub>opt</sub>.160.000 rpm, preferably in a red-contra angle or alternatively in a turbine.

Instruments that are subject to a speed limit of O<sub>max</sub> 300.000 rpm are not authorized for use in a turbine. (e.g. KP6370.314.030/035 or KP6845KR.314.025 – see information on the label)

- Always provide sufficient spray cooling of at least 50 ml/min.
- For good preparation results, the surface should be smoothed with a finisher of matching shape after the initial-preparation.
- The instruments remove substance very effectively. Therefore, only low pressure (<2N) needs to be applied during work.

6 | **DIAO** 

## DIAO. The Range.

#### Instruments for crown preparation Instruments with 8 mm working length



#### Instruments with 10 mm working length



#### Instruments with 12 mm working length

Instruments with 12 mm working length			Occlusal/lingual reduction		
Cylinder, round	Torpedo	Tapered, round	OccluShaper	Egg	Bud
• <b>KP6882L</b> .314.014	• KP6879L.314.014	• KP6856XL.314.021	• KP6370.314.030	• KP6379.314.018	•KP6368.314.016
			• KP6370.314.035	• KP6379.314.023	•KP6368.314. 023
				• <b>KP6379E</b> .314.029	

## **DIAO.** The Range.

#### Instruments with a working length of 6 mm



#### Instruments for cavity preparation



#### Komet Denta

Gebr. Brasseler GmbH & Co. KG Trophagener Weg 25 · 32657 Lemgo Postfach 160 · 32631 Lemgo · Germany

#### Verkauf Deutschland:

Telefon +49 (0) 5261 701-700 Telefax +49 (0) 5261 701-289 info@kometdental.de www.kometdental.de

#### Export

Telefon +49 (0) 5261 701-0 Telefax +49 (0) 5261 701-329 export@kometdental.de www.kometdental.com Komet Austria Handelsagentur GmbH Hellbrunner Straße 15 5020 Salzburg · Austria

Telefon +43 (0) 662 829-434 Telefax +43 (0) 662 829-435 info@kometdental.at www.kometdental.at

Komer, Janz Y, Cerabur Y, Shaper Y, JuShaper Y, OptiPost Y, Polybur Y, PrepMarker Y, Procodile QV, R6 ReziFlow and SHAX" are ademarks of Gebr. Brasseler GmbH & Co. KG.

\* and EndoPilot\* are registered trademarks company Schlumbohm GmbH & Co. KG

iome of the products and designations mentioned in he text are trademarked, patented or copyrighted, he absence of a special reference or the sign <sup>®</sup> should not be terpreted as the absence of legal protection.

This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or reprocessed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

Colours and products are subject to alterations. Printing errors excepted.

As at february 2023

www.kometdental.com/products/diao