For younger patients, historical treatments are not always pulp friendly and show many drawbacks. With Biodentine™, you bring to children the better care they deserve for their dental treatments with an easy and fast procedure for you.

Technical Insights

Higher biocompatibility and bioactivity for pediatric indications

- High biocompatibility assessed and evidenced through 15 scientific publications
- Highest amount of calcium and hydroxide ions released upon setting (1)
- Induces thick dentin bridge formation (2) thanks to largest calcium surface concentration compared to similar dental materials (3)
- Shows both osteogenic and angiogenic properties to promote pulp & tissues healing (4)
- Calcium silicate based materials show better clinical & radiographic results than Formocresol (5)
- Allows the tooth to grow normally

Immediate post-op X-ray after Biodentine placement. The open shows the apical closure. apex is clearly visible.



The 9-month follow-up Xrav


Courtesy Prof. L. Martens & Prof. R. Cauwels, UZ Ghent, Belgium

Bulk fill placement for a strong restoration

- Similar mechanical properties as dentin allowing durable bulk fill procedure
- Fast growing mechanical strength allowing solid restoration as soon as Biodentine[™] is set
- Shows no depth of cure limitation thanks to its biosilicate chemistry

Product properties designed for pediatric dentistry

- Short setting time of 12 min allowing rapid placement of a stainless steel crown when necessary
- No post-operative pain (6)
- No tooth discoloration (7)
- Antimicrobial properties thanks to alkaline ph(=12)(8)



Total handling time	
min	
Setting time in mouth	
6 min	

Kurun Aksoy M, Tulga Oz F, Evaluation of calcium (Ca2+) and hydroxide (OH-) ion International Journal of 2017 diffusion rates of indirect pulp capping materials Orhan K. Artificial Organs Nanoscale chemical surface characterization of Gong V, França R. Journal of Dentistry four different types of dental pulp-capping materials Costa F, Sousa Gomes P, Osteogenic and Angiogenic Response to 2016 Journal of Endodontics Calcium Silicate-based Endodontic Sealers Fernandes MH. El Meligy OA, Allazzam S, Comparison between Biodentine[™] and formocresol for Quintessence Alamoudi NM. pulpotomy of primary teeth: A randomized clinical trial Comparative evaluation of calcium silicate- based dentin substitute Grewal N. Salhan R. (Biodentine[™]) and calcium hydroxide (pulpdent) in the formation of Contemporary Clinical Kaur N. Patel HB reactive dentin bridge in regenerative pulpotomy of vital primary Dentistry teeth: Triple blind, randomized clinical trial Özyürek T, Comparison of the antimicrobial activity of direct pulp-capping Journal of Conservative materials: Mineral trioxide aggregate-Angelus and Biodentine™ Demiryürek EÖ. Dentistry Martens L. Endodontic treatment of trauma-induced necrotic immature European Journal of Rajasekharan S, teeth using a tricalcium silicate-based bioactive cement. Pediatric Dentistry A report of 3 cases with 24-month follow-up. Cauwels R. Evren OK, Altunsoy M, Fracture resistance of simulated immature teeth after European Journal of Tanriver M, Capar ID, 2016 apexification with calcium silicate-based materials Pediatric Dentistry Kalkan A, Gok T. Nowicka A, Wilk G, Tomographic Evaluation of Reparative Dentin Formation after Lipski M, Kołecki J, Direct Pulp Capping with Ca(OH), MTA, Biodentine™ Journal of Endodontics 2015 Buczkowska-Radlińska J. and Dentin Bonding System in Human Teeth Staining Potential of Neo MTA Plus, MTA Plus, Camilleri J. Journal of Endodontics 2015 and Biodentine™ Used for Pulpotomy Procedures Koubi G, Colon P, Franquin Clinical evaluation of the performance and safety of a JC, Hartmann A, Richard G, new dentin substitute, Biodentine™, in the restoration of Clinical Oral Investigation 2012 6 Faure MO. Lambert G. posterior teeth - a prospective study Oral Surgery, Oral Peng L1, Ye L, Evaluation of the formocresol versus mineral trioxide aggregate Medicine, Oral Pathology, 2006 5 Tan H, Zhou X primary molar pulpotomy: a meta-analysis Oral Radiology

Presentation

Available in:

- Box of 15 capsules and 15 single-dose containers
- Box of 5 capsules and 5 single-dose containers



Journal

Septodont - 58 Rue du Pont de Créteil - 94100 Saint-Maur-des-Fossés - France Tél.: +33 (0)1 49 76 70 00 - Fax: +33 (0)1 48 85 54 01 Please visit our website for more information:

www.septodont.com









Dentin Substitute



3 Volets Biodentine 0418 Pediatric.indd 1-3

Biodentine[™]: the better standard for younger patients' teeth 5



Maintain pulp vitality

- Biodentine[™] is bioactive and promotes the pulp's self healing
- In indirect and direct pulp caps, in partial or total pulpotomies: Biodentine helps you save the pulp everytime it's not inflamed
- Its bioactivity allows the apical closure of immature teeth
- Its high pH reduces the risk of bacterial proliferation.



Adapted to younger patients

- Biodentine[™] is "bulk filled" into the cavity
- It can be left as a temporary restoration for up to 6 months
- The final restoration in the same session is possible, whether a direct composite or a stainless steel crown
- No post-operative pain & no discoloration
- No need for extra bonding, no extra steps.



Final restoration placed within 6 months

Clinical cases

Pulpotomy

This includes pulpotomy in temporary teeth and partial pulpotomy in permanent teeth. Beside its bioactivity, Biodentine[™] fast setting time allows immediate crown restoration. It also allows to make it directly intraorally functional without fear of material deterioration



Initial clinical view of #55 for a 8-year-



Pulp exposure during the course of caries curettage necessitates carrying out chamber to the occlusal surface. a cervical pulpotomy (vital pulp, non-



Bulk fill Biodentine™ from the pulp



One month later, placement of a paedodontic cap



X-ray follow-up image at 3 months showing the absence of any periradicular lesion.

Traumatized teeth

Incisors are mostly concerned when it comes to crown fractures. Treatment of these complicated fractures in immature teeth is often a real challenge since the goal is to stimulate the pulp so that the apex can pursue its physiological maturity process. Biodentine[™] shows dentin-like mechanical properties, triggers no tooth discoloration and has a short setting time (12min). It is therefore the material of choice for fragilized immature teeth.



Complicated crown fracture with reversible pulpitis of tooth #21.



Partial pulpotomy was carried out, then the hemorrhage was controlled.



Bulk fill Biodentine™ as a pulp capping material and a temporary restoration.



After one week, Biodentine[™] had fully set and had not washed out, the patient was asymptomatic and final restoration was done




image at 12 months showing continued root

Deep caries treatment

After curetage of caries, the pulp may be seen by transparency or may be exposed. Biodentine™ is a restorative material that can directly be placed on the pulp and is evaluated as superior to MTA, CaOH and Formocresol in numerous publications. It offers zero discoloration ot teeth, provides an outstanding seal, is easy to handle, sets fast and is antibacterial. You can therefore "bulk fill" Biodentine™ without extra conditioning and bonding.



deep cavities in close proximity to the pulp. before and after treatment.



7 years old tooth #46 after caries removal, Indirect capping with Biodentine™, Figure A and B show the radiographic image





12-months, showing the continued formation of radicular