


**ruscello**<sup>®</sup>  
Toothbrushes



**./GC./**

# Grappo

# B-30M




Wider brush head helps provide 1.8 times better plaque removal.

### Available in:

- Red
- Purple
- Blue
- Green
- Yellow

Designed for faster and more effective plaque removal, Grappo also has a thick handle and slip-resistant rubber to ensure a better grip.



With a wider brush head and denser duo-level bristles Grappo helps remove plaque quickly and has the ability to brush the upper and lower teeth at the same time.

Duo-level bristles helps reach into the interproximal area with ease.


# OPeration

# OP-10



Designed for patients who have difficulty brushing due to severe pain from periodontitis or post surgery.


With a dome-like design and super fine 4 mil bristles this toothbrush helps minimise irritation to sensitive gums, yet gently and effectively cleans.





A toothbrush dedicated to the maintenance of dental implants.

Made with hybrid duo-level bristles and in a concave design, this toothbrush is able to adapt to the contoured shape of implants, penetrating the space between implants and gums with ease.



A toothbrush with colour coded bristles and a longer handle designed to aid in reinforcing optimal brushing techniques in young children by the dental professional or parent.




Different character designs available.

# Features of ruscello®


## Brush Head

Difference in brush head size



## Neck

The gentle 80mm curve of Ruscello's neck enables it to clean the hard-to-reach last molars, interproximal spaces and distal molar surface with ease.



## Handle

An ergonomically shaped handle that provides a stable, comfortable grip and allows the brush to reach posterior and anterior teeth easily.




## Bristles

Different bristles for different needs.

### Basic

For daily oral care




### Operation

For post surgical operation patients



### Implant

For patients with implants



## Grip

Functional, non-slip, soft rubber grip allows you to hold the toothbrush firmly.



GC Australasia Dental Pty Ltd  
1753 Botany Road Banksmeadow  
NSW 2019 Australia

T: +61 2 9301 8200  
E: info.australasia@gc.dental  
[www.gcaustralasia.com](http://www.gcaustralasia.com)