

TePe Special Care™ Compact

Features

- 9,000 ultra-soft filaments for delicate cleaning.
- First brush stage after oral surgery.
- For patients with extremely sore and delicate tissues as a result of trauma, radiotherapy treatment, or after oral surgery. * * It is always the dental professional who decides when it is suitable to start brushing after oral surgery.
- Neck can be angled for easier reach.
-

Specification

Handle: Polypropylene (ISCC certified PP)/ Filaments: Polyamide (PA)

Packaging material

Product Packaging: Plastic flowpack = Polypropylene (PP)

Recycling

Package: Plastic flowpack: Container for plastic recycling. Please follow local regulations for recycling in your country/municipality.
Product: Household Waste

Packaging language

English



Item number 452151

EAN consumer

EAN Retailer 7317400005751

EAN Transport 7317400005799

EAN Pallet

Amount in pack 1 in consumer pack

Amount in pack 25 in retail pack