

AquaCare

DISCOVER AQUACARE TWIN



WHY THE TWIN?

- ✓ The Difference between just using a device or investing in a Philosophy
- ✓ Double your options and help you in maximising your treatments
- ✓ Minimum effort with its 1 SEC SWITCH for Maximum Clean
- ✓ Versatility: Change powder in a second and adapt to any situation
- ✓ Minimally Invasive

DESCRIPTION

Clean Dentistry. AquaCare via its use of various powders textures the tooth surface in for improved bond strengths, desensitising, remineralising and prophylaxis treatments. Fissures can be cleared, and cavities can be cut – in many cases without the need for the drill or needle. You are in total control of powder, pressure and solution – what treatments might your patients require today with our Supra or Sub Gingival tips...?

SPECIFICATIONS

Adhesive Dentistry, Bonding, Clean, Desensitizing, Endo, Hygiene, Implants, Laboratory, Lithium Disilicate, MI, Ortho, Pediatric, Pekkton®, Perio, Remineralising, Restorative Dentistry, Stain Removal, Veneers, Zirconia Oxide... Supra or Sub Gingivally.

PRODUCT FEATURES

- ✓ 2 Chambers
- ✓ 3 in 1 Handpiece
- ✓ 3 in 1 Foot pedal

ORDERING DETAILS

PART NUMBER: I/MAC8502F
AMOUNT PER ORDER: 1

WARRANTY: 1 Year (extended to 2 after 1 year service)



RELATED PRODUCTS



AquaCare Powder



AquaCare Cart



Dust Cabinet

AquaCare

DISCOVER AQUACARE SINGLE



WHY THE SINGLE?

- ✓ Minimally invasive
- ✓ All the versatility, One procedure at a time
- ✓ One thought in mind: Ideal for the hygienist that wants to Clean with Aluminium Hydroxide, Sodium Bicarbonate, Sodium Bicarbonate Soft or Syc. Or the for the clinician that prefers to sandblast.

DESCRIPTION

Clean Dentistry. AquaCare via its use of various powders textures the tooth surface in for improved bond strengths, desensitising, remineralising and prophylaxis treatments. Fissures can be cleared, and cavities can be cut – in many cases without the need for the drill or needle. You are in total control of powder, pressure and solution – what treatments might your patients require today with our Supra or Sub Gingival tips...? The main differences between the Single and the Twin is time. Time to change the cartridges, time to purge, time to alter the settings.

SPECIFICATIONS

Adhesive Dentistry, Bonding, Clean, Desensitizing, Endo, Hygiene, Implants, Laboratory, Lithium Disilicate, MI, Ortho, Pediatric, Pekkton®, Perio, Remineralising, Restorative Dentistry, Stain Removal, Veneers, Zirconia Oxide... Supra or Sub Gingivally.

PRODUCT FEATURES

- ✓ 1 Chamber
- ✓ 3 in 1 Handpiece
- ✓ 3 in 1 Foot pedal

RELATED PRODUCTS



AquaCare Powder



AquaCare Cart



AquaSol

ORDERING DETAILS

PART NUMBER: I/MAC8502F
AMOUNT PER ORDER: 1

WARRANTY: 1 Year (extended to 2 after 1 year service)

AquaCare

DISCOVER AQUACARE TIPS



Standard Tip

iTip

Plus Tip

WHY OUR RANGE OF TIPS?

- ✓ Precision
- ✓ Excellence
- ✓ Innovation
- ✓ Detail-oriented

DESCRIPTION

Standard Tip: Suitable for all powders and both handpieces.

iTip: Designed for interproximal and diastema use, though routinely used in many other applications **ONLY TO BE USED WITH THE 0.6mm HANDPIECE.**

Plus Tip: Ideal for larger particle size and wider coverage, e.g. Hypomineralised cases & hygiene use.

FEATURES

With more tips in development, each of our tips is adapting to procedures and treatment and helping you to reach 99% precision.

RELATED PRODUCTS



AquaCare Powder



AquaCare Handpieces



AquaCare Twin

ORDERING DETAILS

PART NUMBERS:

Standard Tip: I/ASS8007F
Plus iTip: I/ASS8009F
Plus Tip: I/ASS8008F

AMOUNT PER ORDER: Pack of 50 per type

AquaCare

DISCOVER AQUACARE DRYAIR



WHY AQUACARE DRYAIR?

- ✓ Specifically designed for AquaCare
- ✓ Perfectly DryAir supply for a perfect work with Air-Abrasion and Polishing
- ✓ Reduce the risk of moisture contamination from your air-line

DESCRIPTION

The DryAir removes any moisture contamination from the compressed air-line used to drive the AquaCare unit. An excellent dry air supply is essential for trouble-free operation of AquaCare units.

SPECIFICATIONS

Our compact DryAir air drying unit can be integrated into the stand and will ensure the quality of the air is suitable for air abrasion and air polishing powders.

The DryAir unit is also available as a stand-alone device where a stand is not required.

PRODUCT FEATURES

- ✓ Specifically designed for AquaCare
- ✓ Designed to fit perfectly with you AquaCare Cart

RELATED PRODUCTS



AquaCare Handpieces



AquaCare Cart



AquaCare Twin

ORDERING DETAILS

PART NUMBER: I/MAC8552F
AMOUNT PER ORDER: 1

WARRANTY: 1 Year

AquaCare

DISCOVER AQUACARE DENTAL AIR ABRASION CABINET



WHY THE AQUACARE DENTAL AIR ABRASION CABINET?

- ✓ Perfect for extra-oral work
- ✓ Perfect for lab usage
- ✓ Ensure a 100% clean environment

DESCRIPTION

It offers a simple aid in support of your AquaCare unit to enable abrasion procedures to be done outside of the mouth.

Dental attachments such as crowns, bridges or brackets can be cleaned and abraded before they are cemented in the mouth.

SPECIFICATIONS

The cabinet is very efficient and easy to use. It has been designed to be small and mobile so that it can find a place in the smallest surgery without disturbing the arrangement of your workplace.

It is self-contained and requires no connection to a power supply – it uses a filter foam pad to catch the dust within the cabinet. The abrasive powder and debris retained in the filter can be removed as and when needed and the filter can be reused after cleaning.

Moreover, pieces like the Foam Filter, Side Entry Gloves or the Clear View Window can be replaced as needed when worn.

PRODUCT FEATURES

- ✓ User centred
- ✓ Easy to use and clean
- ✓ Light weight, easy to move around when required

RELATED PRODUCTS



AquaCare Powder



AquaCare Handpieces



AquaCare Twin

ORDERING DETAILS

PART NUMBER: I/MAC6200F
AMOUNT PER ORDER: 1

WARRANTY: 1 Year

AquaCare

DISCOVER AQUACARE POWDERS THE INTRODUCTORY PACK



WHY AQUACARE POWDERS?

- ✓ Superdry Powder
- ✓ Built-in flow control
- ✓ User-friendly with easy colour code cartridges
- ✓ Quality assured
- ✓ Multidisciplinary range

DESCRIPTION

29 µm Al₂O₃ Abrasive powder, for routine dental procedures

53 µm Al₂O₃ Abrasive powder, for hypomineralise and heavier dental procedures

NaHCO₃ Cleaning powder, effective in the removal of the biofilm and general prophy cleaning

Sylc BAG45S5 Bioglass, Desensitising, Cleaning powder Remineralising and desensitising, stain removal for SUPRA & SUB GINGIVAL

SPECIFICATIONS

We've designed AquaCare so that it works ideally and allows for the broadest range of dental cases.

That's why we invented a cartridge system: they're easy to switch, replace and as more powders become available they are future proof.

RELATED PRODUCTS



AquaCare Tips



AquaCare Handpieces



AquaCare Twin

ORDERING DETAILS

PART NUMBERS:

29 µm Aluminium Oxide: I/PDR8025F ; 53 µm Aluminium Oxide: I/PDR8024F

Sylc - Desensitising: I/PDR0034F ; Sodium Bicarbonate: I/PDR8014F

AMOUNT PER ORDER: Box of 4 per type

AquaCare

DISCOVER AQUACARE AQUASOL



WHY AQUACARE AQUASOL?

- ✓ Keeps the mess away
- ✓ Vanilla taste for patient comfort
- ✓ Designed with a fine amount of ethanol to auto-clean your AquaCare while working (prevents the build up of bacteria)

DESCRIPTION

AquaSol Fluid is ideal for use with the AquaCare as the fluid component to AquaAbrasion™. With a beautiful vanilla fragrance, your patients will enjoy their AquaCare treatment. A 17.5% ethanol content helps to sterilise the inner workings of your AquaCare, helping to maintain hygienic procedures.

SPECIFICATIONS

Bottle Volume: 500ml / 16.9floz
Composition: Water based solution containing 17.5% Ethanol
Fragrance: Vanilla
Suitable for: All AquaCare Units

RELATED PRODUCTS



AquaCare Powder



AquaCare Handpieces



AquaCare Twin

ORDERING DETAILS

PART NUMBERS:

Per 1 Bottle: I/PAC8200F
Per Pack of 6 Bottles: I/PAC8206F

AMOUNT PER ORDER: 1 or Box of 6.

AquaCare

DISCOVER AQUACARE HAND PIECES



WHY AQUACARE HAND PIECES?

- ✓ Light weight
- ✓ Precise control of spray
- ✓ Auto-clave
- ✓ Never clog

DESCRIPTION

The AquaCare's patented handpiece delivers a high-speed fluid stream, fed from the AquaSol, to create a liquid curtain around the abrasive medium.

The autoclavable handpiece comes in colour-coded nozzle sizes. The long-life tungsten carbide nozzles are easily changed to suit the treatment.

An entire procedure can be performed without removing the tip from the oral cavity. One handpiece supplies air, water and treating medium.

SPECIFICATIONS

Lightweight and fully autoclavable hp and sheaths available in 0.06mm Silver and 0.08mm Gold

RELATED PRODUCTS



AquaCare Tips



AquaCare Powders



AquaCare Twin

ORDERING DETAILS

PART NUMBERS:

Silver Hand Piece: I/ASS8191F
Gold Hand Piece: I/ASS8192F

AMOUNT PER ORDER: 1 per type