

# Think of it as a routine check up

for your steriliser.

The **ENSURE**<sup>™</sup> line of infection control monitoring products confirms your equipment is performing at its peak. It includes a range of products that are simple, easy to use and help to ensure your patients are being kept safe.



#### Easy

Provide the highest level of sterility assurance with easy-to-read indicators, which clearly identify if an instrument load met the expected parameters for sterilisation.



#### Safe

Designed to help keep your patients safe.

## Capacity

The 10-well incubator allows you to monitor up to 5 sterilisers at one time, providing you with easy-to-read results that show whether sterility has been achieved or not.

Reliable

Meets the most stringent requirements in today's changing regulatory environment.



#### **ENSURE**<sup>™</sup> Incubator

The ENSURE<sup>™</sup> Incubator is a 10-well incubator with an integrated vial crusher and built-in timer for up to five different sterilisers.



#### **Biological Indicators (BI)**

Fast, reliable, self-contained indicators provide results in as little as 24 hours. The only cost effective method to ensure effective performance of your steriliser.

The ENSURE BI's are available in a starter kit system that consist of 25 BI's, the incubator and an easy to fill-in log book. The starter kit is the most cost effective, accurate system for monitoring your steam autoclave directly in your office.



### Bowie Dick Test (BD) for Steam Penetration

The ENSURE BD test is designed for steam vacuum assisted autoclaves to ensure proper steam penetration for a confirmed sterilisation process. SciCan's BD test is a unique small form factor test that is easy to use and stored. The BD test is available as a starter kit that contains the unique test card holder and 31 test cards.

| ENSURE INCUBATOR |  |
|------------------|--|
| Part #           | Description  |
| SCI-BISK-230EU   | Ensure system, Starter Kit<br>(25 Bl's, incubator, log book) |
| SCI-BINC-230     | Ensure incubator (incubator only)                            |
| SCI-BILOG        | Ensure biological monitoring system log book                 |
|                  |  |

| ENSURE BIOLOGICAL INDICATORS (BI) |                                       |  |
|-----------------------------------|---------------------------------------|--|
| Part #                            | Description                           |  |
| SCI-BI25                          | Ensure BI's (25 count<br>refill box)  |  |
| SCI-BI100                         | Ensure BI's (100 count<br>refill box) |  |
|                                   |                                       |  |

| ENSURE BOWIE DICK TEST (BD) |   |
|-----------------------------|---|
| Part #                      | Description   |
| SCI-BDSK134V                | Ensure Bowie Dick Starter<br>Kit, 134°C (Holder + 31 test<br>cards) |
| SCI-BD134V Class 2          | Ensure Bowie Dick Test<br>Cards (30 count)                          |

Manufactured for: SciCan Ltd. 1440 Don Mills Road Toronto ON M3B 3P9

CANADA Phone: 416 445 1600 Toll free: 1-800-667-7733 Fax: 416 445 2727

EU Rep/Distributed by: SciCan GmbH Wangener Strasse 78 88299 Leutkirch GERMANY Tel.: +49 (0)7561 98343 - 0 Fax: +49 (0)7561 98343 - 699

#### Distributed by: SciCan Medtech AG Alpenstrasse 14 6300 Zug

SWITZERLAND Tel.: +41 (0) 41-727 7027 Fax: +41 (0) 41-727 7029



SciCan Your Infection Control Specialist

www.scican.com

