

Essentials Nitrile Examination Gloves

**Enhanced grip
and protection.**



NEW!



Order today at
henryschein.com.au

 **essentials**
everyday supplies. everyday value.

Essentials Nitrile Examination Gloves

- ✓ Microtextured fingertips ensure a secure grip on instruments and materials
- ✓ Beaded cuff prevents rolling back, 0.07mm thickness at fingertips/0.05mm thickness at palm
- ✓ Powder free to reduce the risk of skin and respiratory sensitivities
- ✓ Soft, comfortable ambidextrous design
- ✓ Single-use
- ✓ Non-sterile
- ✓ Not made with natural rubber latex



CLASS 1 MEDICAL DEVICE CONFORMING TO REGULATION
EU MDR 2017/745 AND COMPLY TO EN 455-1/2/3/4
AND ASTM D6319 • AQL 1,5



Essentials Nitrile Examination Gloves (Blue)

SIZE	QUANTITY	ITEM CODE
XS	100/Box	HSE-5724318-C
S	100/Box	HSE-5724319-C
M	100/Box	HSE-5724320-C
L	100/Box	HSE-5724321-C
XL	100/Box	HSE-5724322-C

Supplied as cartons of 10 boxes



Order today at
henryschein.com.au

 **essentials**
everyday supplies. everyday value.