



# HENRY SCHEIN BRAND INFECTION CONTROL

Henry Schein Brand products offer a comprehensive selection of supplies for your practice needs. Every item has met quality and regulatory standards to earn the Henry Schein Seal of Excellence – your 100% guarantee of satisfaction.

 HENRY SCHEIN®

# PERSONAL PROTECTION GLOVES

## ALOE VERA LATEX EXAMINATION GLOVES



### Moisture for your hands


Powder-free latex gloves containing aloe vera gel which moisturises and hydrates the skin. Excellent tactile sensitivity. Green in colour. Available in sizes extra small to extra large.

Box of 100 10 boxes per carton

XS	HS-9002871
Small	HS-9002872
Medium	HS-9002873
Large	HS-9002874
Extra Large	HS-9002875

### Standards

Complies with AS/NZS 4011.1:2014.  
Manufactured to EN 455-1, 2, 3.  
AQL (water tight): 1.5  
ARTG #121432  
Single use, non sterile



Textured palms and fingertips

## CRITERION CL LATEX EXAMINATION GLOVES



### Enhanced Grip

A quality grade glove designed for a comfortable fit. Chlorinated for easily donning and doffing. Fully textured for optimal grip. Powder free latex.

Box of 100 10 boxes per carton

XS	HS-102-5354
Small	HS-102-8296
Medium	HS-102-5355
Large	HS-102-7763

### Standards

Complies with AS/NZS 4011.1:2014.  
AQL (water tight): 1.5  
ARTG #121432  
Single use, non sterile



Full Texture

## NITRILE EXAMINATION GLOVES SCENTED



### Latex free scented

These unique powder free nitrile gloves offer excellent fit and have a pleasant scent of green tea or bubblegum. Their superior elasticity makes them comfortable to wear, and the finger texture ensures secure grip.

Box of 100. 10 boxes per carton

<b>Green Tea Scented</b>	pale green
Small	HS-9883761
Medium	HS-9883762
Large	HS-9883763
<b>Bubble Gum Scented</b>	pink
Small	HS-9883765
Medium	HS-9883766
Large	HS-9883767

### Standards

EN 455-1,2,3,4. ASTM D6319  
AQL (water tight): 1.5  
ARTG #300031  
Single use, non sterile



Textured palms and fingertips

# PERSONAL PROTECTION MASKS



**Masks must cover both the nose and mouth and folded out fully to cover the chin.**

Given the risk of splash from fluids and aerosols produced during most dental procedures, the recommended level for masks is a Level 2.

Below are the key criteria for a Level 2 mask according to AS 4381:2015.

- Fluid resistance: 120 mm Hg, ASTM F1862/F1862M-13
- Bacterial Filtration Efficiency (BFE%):  $\geq 98\%$ , ASTM F2101-14
- Differential Pressure (mm H<sub>2</sub>O/cm<sup>2</sup>):  $< 5.0$ , EN 14683:2014, Appendix C

## HENRY SCHEIN PROCEDURE LEVEL 2 MASK


Comfortable earloop mask with 3ply construction, latex and fiberglass free. Suitable for use in a healthcare environment as it provides effective filtration and fluid resistance (meets requirements of level 2, AS/NZS 4381:2015).

Blue Earloop HS-988-2733

### Standards

AS 4381:2015, Level 2

ARTG # 321807



## MAXIMA MASKS

MAXIMA Masks are available in earloop and tie on. These procedure masks are intended to be worn by healthcare personnel. They offer high bacterial filtration and are fluid resistant to 120mm Hg. The protection provided is of a Level 2 barrier, in accordance with AS 4381:2015.

The triple layer construction consists of fibreglass-free material surrounded by soft, non-irritating, inner and outer facing. Earloops are ultrasonically sealed for comfort.

Blue Earloop HS-9009061

White Earloop HS-9009062

Pink Earloop HS-9009061


Blue Tie-on HS-9009068

### Standards

AS 4381:2015, Level 2

EN 14683:2019, Type IIR

ARTG # 321807



# SURFACE MANAGEMENT CLEANING

Cleaning with a neutral detergent removes bioburden from surfaces and should be the first line of defence when it comes to surface management.

## CleanND WIPES

Henry Schein CleanND Wipes incorporate a high-quality wipe with a special texture to aid cleaning and leave surfaces streak free. Pre-soaked in a pH neutral, biodegradable detergent solution, these wipes are convenient and safe to use in a healthcare setting.

- Thicker more absorbent wipes stay wet and the special texture helps lift biological contaminants and soil from surfaces.
- Detergent solution utilises a biodegradable formula, safe to use.
- Non-corrosive, non staining and streak free.
- Approved for use on medical devices, equipment and hard surfaces in hospitals and healthcare environments

Canister            **HS-572-2628**  
 Refill:             **HS-5722629**  
 Wipes per tub:    180  
 Refill:             180  
 Wipe quality:     Premium  
 ARTG # 329554



## NEUTRAL DETERGENT WIPES

Henry Schein Neutral Detergent Wipes are a surface cleaner approved for use with medical devices. They can be used for pre-cleaning prior to disinfection or sterilisation procedures.

- Neutral pH, non-ionic solution
- Non staining and biodegradable detergent
- Non-corrosive to instruments and handpieces
- Compatible with disinfection and sterilisation programs
- For use on medical devices

Canister            **HS-NWIPE**  
 Refill:             **HS-NWIPESR**  
 Wipes per tub:    220  
 Refill:             220  
 Wipe quality:     Economy  
 ARTG # 186414



# SURFACE MANAGEMENT DISINFECTION

Henry Schein Surface Disinfectant Wipes are a premium quality wipe soaked in a cleaning and broad-spectrum disinfection solution. Developed to reduce cross contamination in healthcare environments, they provide fast acting disinfection including activity against SARS-CoV-2 (Covid 19) in 2 mins.

- 2 in 1 cleaning and disinfection
- Fast acting, broad spectrum, hard surface disinfectant
- Active against SARS-CoV-2 (Covid 19) in 2 mins
- Non-abrasive, non-corrosive
- Lint free, absorbent wipes aid cleaning

Active against *Pseudomonas aeruginosa*, *Escherichia coli*, *Staphylococcus aureus*, Methicillin Resistant *Staphylococcus aureus* (MRSA) *Enterococcus hirae*, *Acinetobacter baumannii*, *Proteus vulgaris*, *Salmonella enterica*, Poliovirus Type 1, Vancomycin Resistant *Enterococcus faecium* (VRE), *Klebsiella pneumoniae*, SARS-CoV-2 (COVID-19), Adenovirus, Herpes Simplex virus Type 1, *Trichophyton mentagrophytes*, *Candida albicans*, *Candida auris*, *Mycobacterium terrae* (for TB)

Canister **HS-5722626**  
Refill: **HS-5722627**

Wipes per tub: 180  
Refill: 180

Wipe quality: Premium

Hospital Grade Disinfectant tested under dirty conditions. TGO104 & EN Norms.  
Not intended for use on medical devices, therapeutic products or on skin.

ARTG # 335416

CLEANS & DISINFECTS



## ALCOHOL WIPES - ISOPROPYL ALCOHOL 70%

Henry Schein Alcohol Wipes are hospital grade disinfectant with rapid bactericidal activity. To be used on pre-cleaned surfaces in healthcare environments. Use either HS Neutral Detergent Wipes or CleanND Wipes prior. Suitable for use on hard non porous surfaces.

- Hospital Grade hard surface disinfectant
- Bactericidal
- Streak free formulation

Canister **HS-5722624**  
Refill: **HS-5722625**

Wipes per tub: 220  
Refill: 220

Wipe quality: Economy

Compliance to AS/NZS 4187 and AS/NZS 4815.

Hospital Grade Disinfectant. Not intended for use on medical devices, therapeutic products or on skin.

DISINFECTS



# BARRIER PRODUCTS SURFACE MANAGEMENT

Barriers should not be used in lieu of surface management.  
Clean in between barrier change with each patient. ARTG # 311803

NEW

## TRIPLEX SYRINGE SLEEVE



254x64mm  
Pack of 500  
HS-9009372

## SENSOR SLEEVES



35x200mm  
Pack of 500  
HS-9880499



40 x 200mm  
Pack of 500  
HS-9880500

## HANDPIECE SLEEVE



38x203mm  
Box of 500  
HS-9009043

## LED CURING LIGHT SLEEVE



Tip 2.6 cm wide / Body  
6.1 cm wide  
Box of 250  
HS-9880090



## BARRIER FILM ADHESIVE



101x152mm  
Roll of 1200  
Blue HS-9009373  
Clear HS-9009370

## TRAY COVER SLEEVE



36x27cm  
Pack of 500  
HS-9009376

## CHAIR BARRIER SLEEVE



70x61cm  
Pack of 225  
HS-9009374

## MAXI-GARD VL BARRIER ENVELOPE



Size 0 Box of 100  
HS-5701558  
Size 1 Box of 100  
HS-5701559  
Size 2 Box of 300  
HS-5701560

## MAXI-GARD HL BARRIER ENVELOPE



Size 0 Box of 100  
HS-5701561  
Size 1 Box of 100  
HS-5701563  
Size 2 Box of 300  
HS-5701564

## MAXI-GARD MC BARRIER ENVELOPE



Size 0 Box of 200  
HS-5701565  
Size 1 Box of 200  
HS-5701566  
Size 2 Box of 200  
HS-5701567

# STERILISATION CHEMICAL INDICATORS & POUCHES

## MAXIMA SELF SEAL STERILISATION POUCH

Class 1 Chemical Indicator

Maxima Self Seal Sterilisation Pouches are for packaging dental instruments prior to sterilisation and storage once sterilised. These pouches have an external and isolated internal indicator for accurate reading. Designed to indicate that instruments have undergone sterilisation.

Single use Box of 200

57x100 mm (2.25x4")	HS-9006153
70x230 mm (2.75x9")	HS-9006154
90x140 mm (3.5x 5.5")	HS-9006155
90x255 mm (3.5x10")	HS-9006156
133x255 mm (5.25x10")	HS-9006157
190x330 mm (7.5x13")	HS-9006158
230x380 mm (9x15")	HS-9006159

### Standards

Conform to the EN 868-5, ISO 11140-1, and ISO 11607-1.  
ARTG # 115374



## HENRY SCHEIN BOWIE & DICK TEST

Class 2 Chemical Indicator Porous Load


Test set for steam penetration in the autoclave. Colour changes to green to confirm correct function. Lead free and latex free.

Bowie & Dick test kit **9008646**  
Includes metal cassette, card and test sheets

Bowie & Dick test refill **9008647**  
1 card and 20 test sheets. Loading cards can be used 20 times only. Replace loading cards after 20 usages.

### Standards

Conforms to ISO 11140-1



## HENRY SCHEIN HELIX TEST KIT

Class 2 Chemical Indicator Non Porous Load (helix)

Helix steam test strips for autoclaves. Latex-free self-adhesive backing for storage on files. Kit includes Holder plus Sheets x 250.

Helix Test kit **HS-900-8649**  
Includes device and 250 test strips

### Standards

Conforms EN 867-5, ISO 11140 - 1, Class 2.



# SUCTION EVACUATION CLEANERS

NEW

## EUROSEPT XTRA EVAC CLEANER

Effectively dissolves organic and mineral deposits from dental suction lines, keeping them clean and odour free. The EuroSept Xtra Evac Cleaner is a two-component system (Daily and Weekly) for cleaning and maintaining dental suction systems.

To reduce their environmental impact, these products are formulated to be aldehyde-free, phenol-free, and chlorine-free. Regular usage of the DAILY/WEEKLY system helps remove blockages, by dissolving existing organic and mineral deposits, and preventing the build-up of odour causing biofilm in the suction lines.

### Eurosept Xtra Evac Cleaner Daily Concentrate

A low-foaming enzymatic cleaning concentrate, for daily cleaning and care of dental suction systems. Enzymatic cleaners effectively remove organic matter such as blood and saliva residue.

### Eurosept Xtra Evac Cleaner Weekly Concentrate

A foam-free, low pH acid cleaning concentrate, for removal of mineral deposits from dental suction systems. Phosphonates, widely in water treatment processes, is key active in the Weekly Concentrate.

Use DAILY or WEEKLY concentrate at the end of the working day. Note these are to be used on separate days.

The Recommended dilution rate is 1 part in 20 into water (i.e. 5%). The same ratio is used with both DAILY and WEEKLY concentrates.

You will need to combine 100 mL of concentrate and 1.9L of water, to give a total of 2L of the ready to use solution.

### The Range Comprises of:

Eurosept Xtra Evac Cleaner Daily 1L	HS-9888965
Eurosept Xtra Evac Cleaner Daily 5L	HS-9888966
Eurosept Xtra Evac Cleaner Weekly 1L	HS-9888967
ARTG # 384110	



## EUROSEPT XTRA EVAC ACCESSORIES

Eurosept Xtra Dosing Pump for 5L Cannister	HS-9792557
Eurosept Xtra Evac Delivery Bottle	HS-EJUG2L



## SUCTION LINE MANAGEMENT ESSENTIALS WEBINAR

Presented by Emeritus Professor Laurence J. Walsh AO

0.5 CPD Credits

Available on Demand

