

Your **gelatin-based sponge** to control bleeding



Hemogelatin

**Biocompatible & resorbable
hemostatic of porcine origin**



Hemogelatin - gelatin-based hemostatic for efficient bleeding management.

Hemogelatin can be used dry or soaked with sterile saline solution through simple application or by exerting slight pressure like a haemostatic device for controlling the bleeding.

Features	Advantages	Benefits
Gelatin (porcine origin)	Good hemostatic activity ⁽¹⁾ (absorbs 35 times its weight)	Efficient for bleeding management
Resorbable	No need for a second appointment	Optimal patients' flow
Flexible cube shape (10x10x10)	Convenient to use and cut	Easy to place; fits all clinical situations
Individually blister-packed sterile sponges	Maintain the sterility of the product	Peace of mind

Indications

Hemogelatin can be used on:

- Extraction sites
- Periodontal surgical wounds
- Suture sites
- Oral ulcers (non-infected or viral)
- Traumatic wounds
- Denture sores

Product Information

- Porcine origin
- Class III medical device
- Box containing 24 sponges (12 blisters of 2 Individually blister-packed sterile sponges)



(1) internal source file

Distributor:

Septodont - 58 rue du Pont de Créteil - 94107 Saint-Maur-des-Fossés Cedex - France
Tél. : +33 (0)1 49 76 70 00

Please visit our website for more information:

www.septodont.com

