

Articulating Paper

**Simplicity and safety
you can trust.**



Order today at
Henryschein.com.au




EXCLUSIVE BRANDS
AND PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS



Articulating Paper

Henry Schein Articulating Paper

Articulating solutions across booklet, dispenser, refill, and sheet formats

Henry Schein Articulating Paper is intended for identifying occlusal contact points during articulation checks. The range includes booklet, dispenser, refill, and sheet formats across multiple thicknesses to support different clinical preferences, adjustment requirements, and workflow setups.

Product Highlights

- Available in multiple thicknesses, including 38 micron, 40 micron, 45 micron, 71 micron, and 195 micron options
- Formats include booklets, dispensers, refills, sheets, and horseshoe sheets
- Range designed to support different articulation marking preferences and workflow setups
- Suitable for restorative adjustment, occlusal assessment, and general articulation checks

Applications

- Identification of occlusal contact points
- Articulation checks during restorative procedures
- Selection of different thicknesses for varying marking requirements
- Choice of booklet, dispenser, refill, or sheet formats to suit practice workflows

Henry Schein Articulating Paper Range

Choose the right articulating solution for your workflow

Booklets

Product Code	Thickness	Description	Pack Size / Qty
HS-100-8575	38 micron	Articulating Paper Thin Blue	12 Books
HS-100-2876	71 micron	Articulating Paper Thin Blue	12 Books

Dispensers

Product Code	Thickness	Description	Pack Size / Qty
HS-9000725	40 micron	Articulating Paper Dispenser Blue	600 Sheets
HS-9000730	195 micron	Articulating Paper Dispenser Blue	300 Sheets

Refills

Product Code	Thickness	Description	Pack Size / Qty
HS-9000726	40 micron	Articulating Paper Refill Blue	600 Sheets
HS-9000731	195 micron	Articulating Paper Refill Blue	300 Sheets

Sheets

Product Code	Thickness	Description	Pack Size / Qty
HS-9000734	195 micron	Articulating Paper Thick Blue	12 × 12 Sheets
HS-9000735	45 micron	Articulating Paper Blue/Red	6 × 12 Sheets
HS-9000736	45 micron	Articulating Paper Red/Blue	12 × 12 Sheets