

Renfert

making work easy




EASY blank wax

NEW CAD/CAM wax blank

EASY blank wax

CAD/CAM wax blank

CAD/CAM wax blank for the production of digitally designed crowns and bridges for press and casting techniques.



Technical Data

	No. 7700014	No. 7700020
Dimensions (Ø x H)	Ø 98.5 mm x 14 mm (with 10 mm groove)	Ø 98.5 mm x 20 mm (with 10 mm groove)
Dropping point	100–120°C	100–120°C
Weight	approx. 99 g	approx. 138 g
Color	blue	blue

Ordering information

1 EASY blank wax, Ø 98.5 x 14 mm, 1 piece	No. 7700014
2 EASY blank wax, Ø 98.5 x 20 mm, 1 piece	No. 7700020

Advantages

- Precisely milled results with extremely smooth surfaces and strong edge stability due to the high degree of material homogeneity.
- Especially suitable for full ceramic restorations as the wax will burn-out without residue.
- Reliable milling process due to small milled chips which do not stick to the cutter or the object surface.

Making work easy

The Renfert EASY blank wax is the optimal CAD/CAM wax blank for milled restorations in the customarily Renfert quality "made in Germany".

„Making work easy“ stands for outstanding product characteristics, e.g. high material homogeneity. This allows, amongst other things, a perfect fit and surface quality to be achieved in the final restoration.

