

AccuFab-Aris Dental 3D Printer

Just Print It

From Chairside to
Delivery in Minutes.

Empower your clinic with the speed of in-house 3D printing. Start with AccuFab-Aris—where ultimate ease of use enables you to deliver exceptional treatment on the spot. Faster turnover, fewer visits, and happier patients—all from one powerful device.



Compact Design

A Footprint Smaller than Your Laptop.



One-Click Workflow

Precision Printing, Effortlessly Delivered.



Hassle-Free Maintenance

Designed for Non-stop Clinical Productivity.



Solid Performance

Consistent Accuracy for Predictable Results.



Open Material System

Optimized for Shining Dent, Open to Third-party Materials.




Provisional Prosthetic



Denture



Surgical Guide




Splint



Model

A Unified Workflow. Engineered for Clinicians.



Weight	<10KG
Printer Size	240 × 220 x 340 mm
Printer Volume	134 x 84 x 85 mm
Pixel Size	35 μm
Print Speed	Max 80 mm/h
Layer Thickness	50 - 125 μm (Adjustable)
Connectivity	USB / Wi-Fi / Ethernet
Platform Material	Anodized Aluminum
Resin Tank Volume	200 mL
Platform Heating	✓
Resin Tank (NFC) Detection	✓
Residue Detection	✓

*Notice: SHINING 3D reserves the right to introduce modifications or alterations to any of the specifications and images used in this document.

