

# Planmeca Chair™



Planmeca Chair™ provides excellent working ergonomics for the dental team and optimal comfort for patients of all sizes.

## Intelligent design details

*Ergonomic in design, Planmeca Chair™ offers ultimate comfort for both the patient and the dental team.*



### Perfect patient comfort

Planmeca Chair™ is available either with a fixed or automatic leg rest. The *fixed* legrest provides optimal knee angulation and ensures uniform upholstery. Patient entry and exit can be facilitated with the *automatic* foldable leg rest. The perfectly designed backrest and firm armrests maximise patient comfort.

### Easily adjustable headrest

The double articulated headrest allows comfortable positioning of the patient's head providing also good visibility to the treatment area. Thanks to the smooth and lightweight headrest length adjustment, the chair can be easily adapted for patients of all sizes.

### Excellent working ergonomics

The slim and narrow backrest and extensive legroom ensure optimal working ergonomics for the whole dental team. The automatic leg rest facilitates closer face-to-face positioning that is greatly appreciated in doctor-patient consultations and in prosthodontic treatments.



## Options for a variety of needs

*The highly functional Planmeca Chair adapts perfectly to a variety of specialised needs.*



### Flexible positioning

The 180-degree swivel allows the chair to be swiveled 90 degrees left and right. This is a valuable feature when using an intraoral X-ray unit or other auxiliary equipment. The swivel function provides a unique freedom for clinic design and is especially appreciated when space is limited. The mobile base with wheels provides extra possibilities for positioning your Planmeca Chair.

### Great adaptability

The wide movement range enables treating patients freely in various positions. The chair can be adjusted to a very low position, a useful feature when the patient is in a sitting or semi-reclined position. The chair can also be adjusted to a high position level if the dentist prefers to perform a treatment in a standing position.



### Fit for surgical needs

A surgical intravenous armrest is available for immobilising the patient's arm for example during sedation. The surgical headrest option provides enhanced sideways support for the patient's head.



## Unbeatable combination for specialist treatments

You can equip your **Planmeca Chair™** with the unbeatable combination of our large over-the-patient tray and the **Planmeca SingLED®** operating light. This makes the chair an ideal choice for both orthodontic and surgical treatments.

The large and solid over-the-patient tray is made from cast aluminium and has an outstanding load bearing capacity.

The Planmeca SingLED operating light provides optimized illumination to the treatment area and can be adjusted from the upper to the lower jaw position easily.

## High-quality upholsteries for enduring comfort

Planmeca Chair is available with both **Comfy™** and **Ultra Relax™** upholsteries in multiple stylish colours.



# PLANMECA

Asentajankatu 6 | 00880 Helsinki | Finland | tel. +358 20 7795 500 | fax +358 20 7795 555 | sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations. Some products displayed above may not be available in all countries or areas. Rights for changes reserved.

Planmeca, All in one, Anatomat Plus, Cobra, Comfy, Digital perfection, Economat Plus, Elegant, Flexy, Perio Fresh, PlanEasyMill, Planmeca AINO, Planmeca ARA, Planmeca Chair, Planmeca Compact, Planmeca Intra, Planmeca iRomexis, Planmeca Lumion, Planmeca Minea, Planmeca Minendo, Planmeca Minetto, Planmeca Noma, Planmeca Online, Planmeca PlanCAD, Planmeca PlanMill, Planmeca PlanScan, Planmeca Planosil, Planmeca ProCeph, Planmeca ProFace, Planmeca ProID, Planmeca ProMax, Planmeca ProModel, Planmeca ProOne, Planmeca ProScanner, Planmeca ProSensor, Planmeca ProX, Planmeca Romexis, Planmeca SingLED, Planmeca Sovereign, Planmeca Ultra Low Dose, Planmeca Vision, Planmeca Waterline Cleaning System, Proline Dental Stool, Saddle Stool, SmartPan, SmartTouch, Trendy, Ultra Relax and UltraTouch are registered or non-registered trademarks of Planmeca in various countries.

